


POWERHEART® G5

ZOLL®




**DESIGNED FOR
UNEXPECTED HEROES™**



ARE YOU PREPARED TO HELP IN A HEARTBEAT?

It happens in a split second: A person collapses – the victim of sudden cardiac arrest – and the clock starts ticking. You call for help and medical personnel are on their way. But precious seconds are racing by. The victim's best chance for survival depends on immediate CPR and a rescuer with an automated external defibrillator (AED).



CPR plus early defibrillation can **more than double the rate of survival** from out-of-hospital cardiac arrest.

-2015 AHA Guidelines²

THE POWER TO HELP

The latest guidelines from resuscitation councils worldwide are clear: Successful defibrillation must be supported by high-quality CPR.

The AED's first heart analysis is critical. If it calls for no shock, only high-quality CPR can lead to a shockable rhythm on the next heart analysis. High-quality CPR provides a struggling heart with oxygenated blood needed to help restore a normal rhythm.

“

“...monitoring of CPR quality is arguably one of the most significant advances in resuscitation practice in the past 20 years and should be incorporated into every resuscitation and every professional rescuer program.”

- American Heart Association
CPR Quality Consensus Statement³

”

DESIGNED FOR UNEXPECTED HEROES™

The survival rate from sudden cardiac arrest (SCA) is less than 1% worldwide.⁴ According to the American Heart Association (AHA), high-quality CPR, along with early defibrillation, can help improve SCA survival rates.⁵

The Powerheart® G5 AED is the only AED to combine real-time CPR feedback with dual-language functionality, user-paced rescue prompting, and fast shock times.



“

“I was not expecting the AED to walk me through doing compressions. I was sure thankful for that ‘press, press, press’.”

Mark Tucker, Rescuer
CM Tucker Lumber in Pageland, SC

”



Real-Time CPR Feedback

The Powerheart G5 AED uses Intellisense CPR (ICPR) technology to monitor CPR chest compressions, providing rescuers with four distinct corrective voice and text prompts: Press Slower, Press Faster, Press Softer, and Press Harder and Fully Release.



Rescue Ready®

The Powerheart G5 AED performs comprehensive daily self-testing to ensure it is ready to rescue. Studies have shown that approximately 46% of AED failures are due to pad, pad connector, or battery power problems.⁶ The Powerheart AED self-tests detect and help prevent these issues.



User-Paced Rescue Prompts

RescueCoach® prompts guide users step-by-step through the rescue process. The AED recognizes actions taken and ensures that the rescuer completes every critical task before moving on to the next.



Dual-Language Capabilities

With a touch of a button, CPR Instruction, ICPR corrections, and AED prompts can quickly switch from the primary language to a pre-programmed alternate language without delaying the rescue.






THE CASE FOR AEDS

Survival increases with early intervention

Research shows that the probability of survival increases dramatically when CPR is performed and when an AED is applied before an ambulance arrives.¹



Where is the AED?

Too often, the answer is, "We don't have one." Sadly, the same research that demonstrated a nearly fivefold increase in survivability (from 5% to 24%) when an AED is used, also showed that an AED is available only 2% of the time.

POWERHEART® G5



Intellisense™ Defibrillation Pads

Non-polarized pads can be placed in either location on the chest for simple, fast placement.

RescueCoach™

User-paced voice and text prompts guide rescuers through each critical step of the rescue, including CPR instruction.

USB Data Download

90 minutes of rescue data available via USB download.

Dual Language

Change from primary to secondary language with the push of a button, without delaying or interrupting the rescue.



Automatic Self-Testing

Conducts daily self-tests on pads, battery and internal circuitry to ensure your AED is always Rescue Ready®.

Intellisense™ CPR Feedback

Four distinct text and voice corrective CPR prompts for rate, depth and recoil.

Whether you are an experienced rescuer or a first-time responder, you will appreciate how the Powerheart G5 AED combines ease of use, reliability, and advanced technologies to help victims survive sudden cardiac arrest.

The Powerheart G5 AED with Intellisense CPR provides real-time, Guideline-driven feedback to ensure delivery of high-quality compressions. One-button dual-language functionality and daily self-tests give rescuers the confidence to act.

ZOLL MEDICAL CORPORATION

An Asahi Kasei Group Company | 269 Mill Road | Chelmsford, MA 01824 | 978-421-9655 | 800-804-4356 | www.zoll.com
 Copyright © 2020 ZOLL Medical Corporation. All rights reserved. Designed for Unexpected Heroes, Rescue Ready Services, ZOLL, Powerheart, Intellisense and RescueCoach are trademarks or registered trademarks of ZOLL Medical Corporation in the United States and/or other countries. All other trademarks are the property of their respective owners.

¹Weisfeldt ML, et al. J Am Coll Cardiol. 2010;55(16):1713-20.

²AHA 2015 Guidelines: pg. S402

³AHA CPR Quality Consensus Statement. Circulation. 2013; 128; pg. 420

⁴Mehra R; J Electrocardiol. 2007 Nov-Dec;40(6 Suppl):S118-22

⁵Kleinman et al. Part 5: adult basic life support and cardiopulmonary resuscitation quality: 2015 American Heart Association Guidelines Update for Cardiopulmonary Resuscitation and Emergency Cardiovascular Care. Circulation. 2015;132(suppl 2):S414-S435

⁶Annals Emergency Medicine: Analysis of Automated External Defibrillator Device Failures, Reported to the Food and Drug Administration, Volume 59, no. 2, February 2012

For subsidiary addresses and fax numbers, as well as other global locations, please go to www.zoll.com/contacts.

